

CYP2C9 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50815-1**

Capture Antibody Information

Catalog Number: 60576-1-PBS	Clone ID: 1D7G1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cytochrome P450, family 2, subfamily C, polypeptide 9
Isotype: IgG1	GenBank: BC020754	Gene ID: 1559
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag9566	

Detection Antibody Information

Catalog Number: 60576-2-PBS	Clone ID: 1F6G6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cytochrome P450, family 2, subfamily C, polypeptide 9
Isotype: IgG2a	GenBank: BC020754	Gene ID: 1559
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9566	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50815-1 targets CYP2C9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CYP2C9 Monoclonal antibody, PBS Only (Capture) 60576-1-PBS (1D7G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: CYP2C9 Monoclonal antibody, PBS Only (Detector) 60576-2-PBS (1F6G6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50815-1, CYP2C9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60576-1-PBS. Detection antibody: 60576-2-PBS. Standard:Ag9566. Range: 0.098-100 ng/mL.