

REEP4 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50822-1**

Capture Antibody Information

Catalog Number: 60580-1-PBS	Clone ID: 2D2A7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: receptor accessory protein 4
Isotype: IgG1	GenBank: BC013048	Gene ID: 80346
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33149	

Detection Antibody Information

Catalog Number: 60580-2-PBS	Clone ID: 2D4C11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: receptor accessory protein 4
Isotype: IgG1	GenBank: BC013048	Gene ID: 80346
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33149	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP50822-1 targets REEP4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: REEP4 Monoclonal antibody, PBS Only (Capture) 60580-1-PBS (2D2A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: REEP4 Monoclonal antibody, PBS Only (Detector) 60580-2-PBS (2D4C11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50822-1, REEP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60580-1-PBS. Detection antibody: 60580-2-PBS. Standard:Ag33149. Range: 1.563-100 ng/mL.