

RNF133 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50868-2**

Capture Antibody Information

| | | |
|---|---|--|
| Catalog Number: 60611-3-PBS | Clone ID: 1D7D4 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: ring finger protein 133 |
| Isotype: IgG1 | GenBank: BC022038 | Gene ID: 168433 |
| Purification Method: Protein G purification | Immunogen Catalog Number: Ag14534 | |

Detection Antibody Information

| | | |
|--|---|--|
| Catalog Number: 60611-4-PBS | Clone ID: 1A1G12 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: ring finger protein 133 |
| Isotype: IgG3 | GenBank: BC022038 | Gene ID: 168433 |
| Purification Method: Protein A Magarose purification | Immunogen Catalog Number: Ag14534 | |

Applications

| | | |
|--|--|---|
| Tested Applications: Cytometric bead array | Range: 0.098-100 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|--|---|

Product Information

MP50868-2 targets RNF133 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RNF133 Monoclonal antibody, PBS Only (Capture) 60611-3-PBS (1D7D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: RNF133 Monoclonal antibody, PBS Only (Detector) 60611-4-PBS (1A1G12). 100 µg. Concentration 1 mg/mL.

Alternative RNF133 matched antibody pairs: MP00458-1, MP00458-2, MP50868-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50868-2, RNF133 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60611-3-PBS. Detection antibody: 60611-4-PBS. Standard:Ag14534. Range: 0.098-100 ng/mL.