

SEC11A Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50880-3**

Capture Antibody Information

Catalog Number: 67379-5-PBS	Clone ID: 1F9D3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SEC11 homolog A (S. cerevisiae)
Isotype: IgG1	GenBank: BC000359	Gene ID: 23478
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6610	

Detection Antibody Information

Catalog Number: 67379-6-PBS	Clone ID: 2E1E2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SEC11 homolog A (S. cerevisiae)
Isotype: IgG1	GenBank: BC000359	Gene ID: 23478
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6610	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-12.5 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP50880-3 targets SEC11A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEC11A Monoclonal antibody, PBS Only (Capture) 67379-5-PBS (1F9D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: SEC11A Monoclonal antibody, PBS Only (Detector) 67379-6-PBS (2E1E2). 100 µg. Concentration 1 mg/mL.

Alternative SEC11A matched antibody pairs: MP50880-1, MP50880-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50880-3, SEC11A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67379-5-PBS. Detection antibody: 67379-6-PBS. Standard:null. Range: null.

Cytometric bead array standard curve of MP50880-3, SEC11A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67379-5-PBS. Detection antibody: 67379-6-PBS. Standard:Ag6610. Range: 0.098-12.5 ng/mL