

# ZBTB44 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50884-3**

## Capture Antibody Information

<b>Catalog Number:</b> 60621-4-PBS	<b>Clone ID:</b> 2B3F10	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> zinc finger and BTB domain containing 44
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC071729	<b>Gene ID:</b> 29068
<b>Purification Method:</b> Protein G purification	<b>Immunogen Catalog Number:</b> Ag22369	

## Detection Antibody Information

<b>Catalog Number:</b> 60621-5-PBS	<b>Clone ID:</b> 2C3C11	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> zinc finger and BTB domain containing 44
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC071729	<b>Gene ID:</b> 29068
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag22369	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 1.563-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50884-3 targets ZBTB44 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZBTB44 Monoclonal antibody, PBS Only (Capture) 60621-4-PBS (2B3F10). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZBTB44 Monoclonal antibody, PBS Only (Detector) 60621-5-PBS (2C3C11). 100 µg. Concentration 1 mg/mL.

Alternative ZBTB44 matched antibody pairs: MP50884-1, MP50884-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50884-3, ZBTB44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60621-4-PBS. Detection antibody: 60621-5-PBS. Standard:Ag22369. Range: 1.563-100 ng/mL.