

# RABL3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50900-2**

## Capture Antibody Information

**Catalog Number:**  
68342-2-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2G2G2

**Reactivity:**  
human

**GenBank:**  
BC020832

**Immunogen Catalog Number:**  
Ag33137

**Conjugate:**  
Unconjugated

**Full name:**  
RAB, member of RAS oncogene family-Like 3

**Gene ID:**  
285282

## Detection Antibody Information

**Catalog Number:**  
68342-4-PBS

**Host:**  
Mouse

**Isotype:**  
IgG2b

**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
2F4E10

**Reactivity:**  
human

**GenBank:**  
BC020832

**Immunogen Catalog Number:**  
Ag33137

**Conjugate:**  
Unconjugated

**Full name:**  
RAB, member of RAS oncogene family-Like 3

**Gene ID:**  
285282

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50900-2 targets RABL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RABL3 Monoclonal antibody, PBS Only (Capture) 68342-2-PBS (2G2G2). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: RABL3 Monoclonal antibody, PBS Only (Detector) 68342-4-PBS (2F4E10). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50900-2, RABL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68342-2-PBS. Detection antibody: 68342-4-PBS. Standard:Ag33137. Range: 0.098-100 ng/mL.