

RAGE Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50901-1**

Capture Antibody Information

Catalog Number: 67437-1-PBS	Clone ID: 2C8D6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: renal tumor antigen
Isotype: IgG1	GenBank: BC053536	Gene ID: 5891
Purification Method: Protein G purification	Immunogen Catalog Number: Ag29406	

Detection Antibody Information

Catalog Number: 67437-2-PBS	Clone ID: 1C9D11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: renal tumor antigen
Isotype: IgG1	GenBank: BC053536	Gene ID: 5891
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29406	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50901-1 targets RAGE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAGE Monoclonal antibody, PBS Only (Capture) 67437-1-PBS (2C8D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RAGE Monoclonal antibody, PBS Only (Detector) 67437-2-PBS (1C9D11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50901-1, RAGE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67437-1-PBS. Detection antibody: 67437-2-PBS. Standard:Ag29406. Range: 0.781-50 ng/mL