

# NFE2L1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50921-1**

## Capture Antibody Information

**Catalog Number:**  
68634-2-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G purification

**Clone ID:**  
1C4E3

**Reactivity:**  
human

**GenBank:**  
BC010623

**Immunogen Catalog Number:**  
Ag32992

**Conjugate:**  
Unconjugated

**Full name:**  
nuclear factor (erythroid-derived 2)-like 1

**Gene ID:**  
4779

## Detection Antibody Information

**Catalog Number:**  
68634-3-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2F11A6

**Reactivity:**  
human

**GenBank:**  
BC010623

**Immunogen Catalog Number:**  
Ag32992

**Conjugate:**  
Unconjugated

**Full name:**  
nuclear factor (erythroid-derived 2)-like 1

**Gene ID:**  
4779

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50921-1 targets NFE2L1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFE2L1 Monoclonal antibody, PBS Only (Capture) 68634-2-PBS (1C4E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFE2L1 Monoclonal antibody, PBS Only (Detector) 68634-3-PBS (2F11A6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50921-1, NFE2L1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68634-2-PBS. Detection antibody: 68634-3-PBS. Standard:Ag32992. Range: 0.195-50 ng/mL