

KIAA0319L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50926-3**

Capture Antibody Information

| | | |
|--|---|------------------------------------|
| Catalog Number: 68717-2-PBS | Clone ID: 2D10C6 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: KIAA0319-like |
| Isotype: IgG1 | GenBank: BC031672 | Gene ID: 79932 |
| Purification Method: Protein G Magarose purification | Immunogen Catalog Number: Ag15263 | |

Detection Antibody Information

| | | |
|--|---|------------------------------------|
| Catalog Number: 68717-6-PBS | Clone ID: 3D5A2 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: KIAA0319-like |
| Isotype: IgG1 | GenBank: BC031672 | Gene ID: 79932 |
| Purification Method: Protein G Magarose purification | Immunogen Catalog Number: Ag15263 | |

Applications

| | | |
|--|--|---|
| Tested Applications: Cytometric bead array | Range: 0.781-100 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|--|---|

Product Information

MP50926-3 targets KIAA0319L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA0319L Monoclonal antibody, PBS Only (Capture) 68717-2-PBS (2D10C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: KIAA0319L Monoclonal antibody, PBS Only (Detector) 68717-6-PBS (3D5A2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50926-3, KIAA0319L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68717-2-PBS. Detection antibody: 68717-6-PBS. Standard:Ag15263. Range: 0.781-100 ng/mL