

SPAM1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50939-2**

Capture Antibody Information

Catalog Number: 60655-3-PBS	Clone ID: 2B5H11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SPAM1
Isotype: IgG2b	GenBank: BC026163	Gene ID: 6677
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag1682	

Detection Antibody Information

Catalog Number: 60655-2-PBS	Clone ID: 1G1E7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SPAM1
Isotype: IgG1	GenBank: BC026163	Gene ID: 6677
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1682	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50939-2 targets SPAM1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPAM1 Monoclonal antibody, PBS Only (Capture) 60655-3-PBS (2B5H11). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPAM1 Monoclonal antibody, PBS Only (Detector) 60655-2-PBS (1G1E7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50939-2, SPAM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60655-3-PBS. Detection antibody: 60655-2-PBS. Standard:Ag1682. Range: 3.125-100 ng/mL.