

# CCDC21 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50940-2**

## Capture Antibody Information

Catalog Number: 60656-3-PBS	Clone ID: 1D9B8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: coiled-coil domain containing 21
Isotype: IgG1	GenBank: BC064528	Gene ID: 64793
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33073	

## Detection Antibody Information

Catalog Number: 60656-2-PBS	Clone ID: 1H9A7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: coiled-coil domain containing 21
Isotype: IgG2a	GenBank: BC064528	Gene ID: 64793
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag33073	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	---------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP50940-2 targets CCDC21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC21 Monoclonal antibody, PBS Only (Capture) 60656-3-PBS (1D9B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCDC21 Monoclonal antibody, PBS Only (Detector) 60656-2-PBS (1H9A7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50940-2, CCDC21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60656-3-PBS. Detection antibody: 60656-2-PBS. Standard:Ag33073. Range: 0.781-100 ng/mL.