

ELAVL2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50972-1**

Capture Antibody Information

Catalog Number:
67097-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2D8C10
Reactivity:
human, mouse, rat, pig, rabbit
GenBank:
BC030692
Immunogen Catalog Number:
Ag28310

Conjugate:
Unconjugated
Full name:
ELAV (embryonic lethal, abnormal vision, Drosophila)-like 2 (Hu antigen B)
Gene ID:
1993

Detection Antibody Information

Catalog Number:
67097-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
1H1A3
Reactivity:
human
GenBank:
BC030692
Immunogen Catalog Number:
Ag28310

Conjugate:
Unconjugated
Full name:
ELAV (embryonic lethal, abnormal vision, Drosophila)-like 2 (Hu antigen B)
Gene ID:
1993

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50972-1 targets ELAVL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ELAVL2 Monoclonal antibody, PBS Only (Capture) 67097-1-PBS (2D8C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: ELAVL2 Monoclonal antibody, PBS Only (Detector) 67097-2-PBS (1H1A3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50972-1, ELAVL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67097-1-PBS. Detection antibody: 67097-2-PBS. Standard:Ag28310. Range: 0.391-100 ng/mL.