

FKBP10 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50974-1**

Capture Antibody Information

| | | |
|---|---|---|
| Catalog Number: 68085-2-PBS | Clone ID: 1F12A11 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: FK506 binding protein 10, 65 kDa |
| Isotype: IgG1 | GenBank: BC016467 | Gene ID: 60681 |
| Purification Method: Protein G purification | Immunogen Catalog Number: Ag31628 | |

Detection Antibody Information

| | | |
|--|---|---|
| Catalog Number: 68085-3-PBS | Clone ID: 1E7H11 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: FK506 binding protein 10, 65 kDa |
| Isotype: IgG1 | GenBank: BC016467 | Gene ID: 60681 |
| Purification Method: Protein G Magarose purification | Immunogen Catalog Number: Ag31628 | |

Applications

| | | |
|--|---|---|
| Tested Applications: Cytometric bead array | Range: 0.098-25 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|---|---|

Product Information

MP50974-1 targets FKBP10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FKBP10 Monoclonal antibody, PBS Only (Capture) 68085-2-PBS (1F12A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: FKBP10 Monoclonal antibody, PBS Only (Detector) 68085-3-PBS (1E7H11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50974-1, FKBP10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68085-2-PBS. Detection antibody: 68085-3-PBS. Standard:Ag31628. Range: 0.098-25 ng/mL