

AXL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50990-2**

Capture Antibody Information

Catalog Number: 68108-2-PBS	Clone ID: 1A4D11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: AXL receptor tyrosine kinase
Isotype: IgG1	GenBank: BC032229	Gene ID: 558
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag14401	

Detection Antibody Information

Catalog Number: 68108-4-PBS	Clone ID: 1E10D3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: AXL receptor tyrosine kinase
Isotype: IgG1	GenBank: BC032229	Gene ID: 558
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag14401	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50990-2 targets AXL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AXL Monoclonal antibody, PBS Only (Capture) 68108-2-PBS (1A4D11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: AXL Monoclonal antibody, PBS Only (Detector) 68108-4-PBS (1E10D3). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50990-2, AXL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68108-2-PBS. Detection antibody: 68108-4-PBS. Standard:Ag14401. Range: 0.195-100 ng/mL.