

NR4A3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51007-3**

Capture Antibody Information

Catalog Number: 60697-4-PBS	Clone ID: 3D7D9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: nuclear receptor subfamily 4, group A, member 3
Isotype: IgG1	GenBank: NM_006981	Gene ID: 8013
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28956	

Detection Antibody Information

Catalog Number: 60697-2-PBS	Clone ID: 1F2A3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: nuclear receptor subfamily 4, group A, member 3
Isotype: IgG3	GenBank: NM_006981	Gene ID: 8013
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag28956	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP51007-3 targets NR4A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NR4A3 Monoclonal antibody, PBS Only (Capture) 60697-4-PBS (3D7D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: NR4A3 Monoclonal antibody, PBS Only (Detector) 60697-2-PBS (1F2A3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51007-3, NR4A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60697-4-PBS. Detection antibody: 60697-2-PBS. Standard:Ag28956. Range: 0.098-100 ng/mL.