

IGFBP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51029-4**

Capture Antibody Information

Catalog Number: 60716-1-PBS	Clone ID: 1F9B7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: IGF binding protein 1
Isotype: IgG1	GenBank: BC057806	Gene ID: 3484
Purification Method: Protein G purification	Immunogen Catalog Number: Ag5852	

Detection Antibody Information

Catalog Number: 60717-3-PBS	Clone ID: 3C3C11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: IGF binding protein 1
Isotype: IgG1	GenBank: NM_000596.4	Gene ID: 3484
Purification Method: Protein G Magarose purification		

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51029-4 targets IGFBP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IGFBP1 Monoclonal antibody, PBS Only (Capture) 60716-1-PBS (1F9B7). 100 µg. Concentration 1 mg/mL.

Detection antibody: IGFBP1 Monoclonal antibody, PBS Only (Detector) 60717-3-PBS (3C3C11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51029-4, IGF1BP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60716-1-PBS. Detection antibody: 60717-3-PBS. Standard:Eg0951. Range: 0.098-100 ng/mL.