

# LAMB1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51088-1**

## Capture Antibody Information

**Catalog Number:**  
67705-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2C2C7  
**Reactivity:**  
human  
**GenBank:**  
BC113455  
**Immunogen Catalog Number:**  
Ag19182

**Conjugate:**  
Unconjugated  
**Full name:**  
laminin, beta 1  
**Gene ID:**  
3912

## Detection Antibody Information

**Catalog Number:**  
67705-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2C1B4  
**Reactivity:**  
human  
**GenBank:**  
BC113455  
**Immunogen Catalog Number:**  
Ag19182

**Conjugate:**  
Unconjugated  
**Full name:**  
laminin, beta 1  
**Gene ID:**  
3912

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51088-1 targets LAMB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LAMB1 Monoclonal antibody, PBS Only (Capture) 67705-2-PBS (2C2C7). 100 µg. Concentration 1 mg/mL.

Detection antibody: LAMB1 Monoclonal antibody, PBS Only (Detector) 67705-3-PBS (2C1B4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:null. Range: null.

Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:Ag19182. Range: 0.195-100 ng/mL.