

# IGF2BP3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51098-2**

## Capture Antibody Information

**Catalog Number:**  
68470-4-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1B11D8  
**Reactivity:**  
human  
**GenBank:**  
BC065269  
**Immunogen Catalog Number:**  
Ag6349

**Conjugate:**  
Unconjugated  
**Full name:**  
insulin-like growth factor 2 mRNA binding protein 3  
**Gene ID:**  
10643

## Detection Antibody Information

**Catalog Number:**  
68470-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1E12E12  
**Reactivity:**  
human  
**GenBank:**  
BC065269  
**Immunogen Catalog Number:**  
Ag6349

**Conjugate:**  
Unconjugated  
**Full name:**  
insulin-like growth factor 2 mRNA binding protein 3  
**Gene ID:**  
10643

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51098-2 targets IGF2BP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IGF2BP3 Monoclonal antibody, PBS Only (Capture) 68470-4-PBS (1B11D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: IGF2BP3 Monoclonal antibody, PBS Only (Detector) 68470-3-PBS (1E12E12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP51098-2, IGF2BP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68470-4-PBS. Detection antibody: 68470-3-PBS. Standard:Ag6349. Range: 0.781-50 ng/mL