

SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51103-2**

Capture Antibody Information

Catalog Number:
60774-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1E3D2
Reactivity:
human
GenBank:
NM_003070
Immunogen Catalog Number:
Ag33074

Conjugate:
Unconjugated
Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2
Gene ID:
6595

Detection Antibody Information

Catalog Number:
60774-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3H1G8
Reactivity:
human
GenBank:
NM_003070
Immunogen Catalog Number:
Ag33074

Conjugate:
Unconjugated
Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2
Gene ID:
6595

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51103-2 targets SMARCA2-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA2-Specific Monoclonal antibody, PBS Only (Capture) 60774-3-PBS (1E3D2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SMARCA2-Specific Monoclonal antibody, PBS Only (Detector) 60774-2-PBS (3H1G8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51103-2, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-3-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 0.391-100 ng/mL