

RAB26 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51188-1**

Capture Antibody Information

Catalog Number: 60816-1-PBS	Clone ID: 1G11B7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: RAB26, member RAS oncogene family
Isotype: IgG1	GenBank: BC066913	Gene ID: 25837
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6463	

Detection Antibody Information

Catalog Number: 60816-2-PBS	Clone ID: 2C1G4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: RAB26, member RAS oncogene family
Isotype: IgG1	GenBank: BC066913	Gene ID: 25837
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6463	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP51188-1 targets RAB26 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAB26 Monoclonal antibody, PBS Only (Capture/Detector) 60816-1-PBS (1G11B7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RAB26 Monoclonal antibody, PBS Only (Detector) 60816-2-PBS (2C1G4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51188-1, RAB26 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60816-1-PBS. Detection antibody: 60816-2-PBS. Standard:Ag6463. Range: 0.195-25 ng/mL