

RBL2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51189-1**

Capture Antibody Information

Catalog Number: 60817-1-PBS	Clone ID: 3D5C12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: retinoblastoma-like 2 (p130)
Isotype: IgG1	GenBank: BC034490	Gene ID: 5934
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag3377	

Detection Antibody Information

Catalog Number: 60817-2-PBS	Clone ID: 3G11A8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: retinoblastoma-like 2 (p130)
Isotype: IgG1	GenBank: BC034490	Gene ID: 5934
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag3377	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51189-1 targets RBL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBL2 Monoclonal antibody, PBS Only (Capture) 60817-1-PBS (3D5C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBL2 Monoclonal antibody, PBS Only (Detector) 60817-2-PBS (3G11A8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51189-1, RBL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60817-1-PBS. Detection antibody: 60817-2-PBS. Standard:Ag3377. Range: 1.563-100 ng/mL.