

For Research Use Only

# Midkine Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51205-1

## Capture Antibody Information

Catalog Number:  
60759-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1G8G2

Reactivity:  
human

GenBank:  
BC011704

Immunogen Catalog Number:  
Ag33793

Conjugate:  
Unconjugated

Full name:  
midkine (neurite growth-promoting factor 2)

Gene ID:  
4192

## Detection Antibody Information

Catalog Number:  
60759-2-PBS

Host:  
Mouse

Isotype:  
IgG3

Purification Method:  
Protein A Magarose purification

Clone ID:  
2A7F2

Reactivity:  
human

GenBank:  
BC011704

Immunogen Catalog Number:  
Ag33793

Conjugate:  
Unconjugated

Full name:  
midkine (neurite growth-promoting factor 2)

Gene ID:  
4192

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51205-1 targets Midkine in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Midkine Monoclonal antibody, PBS Only (Capture) 60759-1-PBS (1G8G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Midkine Monoclonal antibody, PBS Only (Detector) 60759-2-PBS (2A7F2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51205-1, Midkine Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60759-1-PBS. Detection antibody: 60759-2-PBS. Standard:Ag33793. Range: 0.781-100 ng/mL.