

VWA1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51222-2**

Capture Antibody Information

Catalog Number: 60826-1-PBS	Clone ID: 2D7C11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: von Willebrand factor A domain containing 1
Isotype: IgG1	GenBank: BC059409	Gene ID: 64856
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag37064	

Detection Antibody Information

Catalog Number: 60826-3-PBS	Clone ID: 1B10A11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: von Willebrand factor A domain containing 1
Isotype: IgG2a	GenBank: BC059409	Gene ID: 64856
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag37064	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51222-2 targets VWA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VWA1 Monoclonal antibody, PBS Only (Capture) 60826-1-PBS (2D7C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: VWA1 Monoclonal antibody, PBS Only (Detector) 60826-3-PBS (1B10A11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51222-2, VWA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60826-1-PBS. Detection antibody: 60826-3-PBS. Standard:Ag37064. Range: 0.781-100 ng/mL.