

RPGR Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51225-1**

Capture Antibody Information

Catalog Number:
60830-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
4B7C6
Reactivity:
human
GenBank:
BC031624
Immunogen Catalog Number:
Ag27524

Conjugate:
Unconjugated
Full name:
retinitis pigmentosa GTPase regulator
Gene ID:
6103

Detection Antibody Information

Catalog Number:
60830-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2B11H5
Reactivity:
human
GenBank:
BC031624
Immunogen Catalog Number:
Ag27524

Conjugate:
Unconjugated
Full name:
retinitis pigmentosa GTPase regulator
Gene ID:
6103

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51225-1 targets RPGR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RPGR Monoclonal antibody, PBS Only (Capture) 60830-1-PBS (4B7C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RPGR Monoclonal antibody, PBS Only (Detector) 60830-2-PBS (2B11H5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51225-1, RPGR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60830-1-PBS. Detection antibody: 60830-2-PBS. Standard:Ag27524. Range: 0.195-100 ng/mL.