

# UNC93B1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51235-2**

## Capture Antibody Information

<b>Catalog Number:</b> 60836-3-PBS	<b>Clone ID:</b> 3F10E3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> unc-93 homolog B1 (C. elegans)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC101568	<b>Gene ID:</b> 81622
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag27937	

## Detection Antibody Information

<b>Catalog Number:</b> 60836-4-PBS	<b>Clone ID:</b> 4H1G10	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> unc-93 homolog B1 (C. elegans)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC101568	<b>Gene ID:</b> 81622
<b>Purification Method:</b> Protein G purification	<b>Immunogen Catalog Number:</b> Ag27937	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP51235-2 targets UNC93B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UNC93B1 Monoclonal antibody, PBS Only (Capture) 60836-3-PBS (3F10E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: UNC93B1 Monoclonal antibody, PBS Only (Detector) 60836-4-PBS (4H1G10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP51235-2, UNC93B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60836-3-PBS. Detection antibody: 60836-4-PBS. Standard:Ag27937. Range: 0.391-100 ng/mL