

DLK1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51236-2**

Capture Antibody Information

Catalog Number: 67923-4-PBS	Clone ID: 3B11D5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: delta-like 1 homolog (Drosophila)
Isotype: IgG1	GenBank: BC007741	Gene ID: 8788
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag30916	

Detection Antibody Information

Catalog Number: 67923-3-PBS	Clone ID: 2G10B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: delta-like 1 homolog (Drosophila)
Isotype: IgG1	GenBank: BC007741	Gene ID: 8788
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag30916	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51236-2 targets DLK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DLK1 Monoclonal antibody, PBS Only (Capture) 67923-4-PBS (3B11D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: DLK1 Monoclonal antibody, PBS Only (Detector) 67923-3-PBS (2G10B5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51236-2, DLK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67923-4-PBS. Detection antibody: 67923-3-PBS. Standard:Ag30916. Range: 1.563-200 ng/mL.