

# RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51247-3**

## Capture Antibody Information

Catalog Number: 67502-5-PBS	Clone ID: 2B12A6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Rho-related BTB domain containing 3
Isotype: IgG1	GenBank: BC041337	Gene ID: 22836
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5782	

## Detection Antibody Information

Catalog Number: 67502-6-PBS	Clone ID: 1B3F3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Rho-related BTB domain containing 3
Isotype: IgG1	GenBank: BC041337	Gene ID: 22836
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5782	

## Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP51247-3 targets RHOBTB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RHOBTB3 Monoclonal antibody, PBS Only (Capture) 67502-5-PBS (2B12A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RHOBTB3 Monoclonal antibody, PBS Only (Detector) 67502-6-PBS (1B3F3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP51247-3, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-5-PBS. Detection antibody: 67502-6-PBS. Standard:Ag5782. Range: 3.125-100 ng/mL