

DDO Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51260-1**

Capture Antibody Information

Catalog Number: 60850-1-PBS	Clone ID: 2B8C3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: D-aspartate oxidase
Isotype: IgG2a	GenBank: BC032786	Gene ID: 8528
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29876	

Detection Antibody Information

Catalog Number: 60850-2-PBS	Clone ID: 2G8C8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: D-aspartate oxidase
Isotype: IgG1	GenBank: BC032786	Gene ID: 8528
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29876	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51260-1 targets DDO in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDO Monoclonal antibody, PBS Only (Capture) 60850-1-PBS (2B8C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: DDO Monoclonal antibody, PBS Only (Detector) 60850-2-PBS (2G8C8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51260-1, DDO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60850-1-PBS. Detection antibody: 60850-2-PBS. Standard:Ag29876. Range: 0.781-100 ng/mL.