

CLCN1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51269-1**

Capture Antibody Information

Catalog Number: 60857-1-PBS	Clone ID: 1D9D11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chloride channel 1, skeletal muscle
Isotype: IgG1	GenBank: BC112156	Gene ID: 1180
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag13371	

Detection Antibody Information

Catalog Number: 60857-2-PBS	Clone ID: 1C3H2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chloride channel 1, skeletal muscle
Isotype: IgG1	GenBank: BC112156	Gene ID: 1180
Purification Method: Protein G purification	Immunogen Catalog Number: Ag13371	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51269-1 targets CLCN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CLCN1 Monoclonal antibody, PBS Only (Capture) 60857-1-PBS (1D9D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CLCN1 Monoclonal antibody, PBS Only (Detector) 60857-2-PBS (1C3H2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51269-1, CLCN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60857-1-PBS. Detection antibody: 60857-2-PBS. Standard:Ag13371. Range: 0.781-100 ng/mL.