

PD-ECGF Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51270-1**

Capture Antibody Information

Catalog Number: 67104-2-PBS	Clone ID: 1B3E7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymidine phosphorylase
Isotype: IgG2a	GenBank: BC018160	Gene ID: 1890
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28631	

Detection Antibody Information

Catalog Number: 67104-3-PBS	Clone ID: 1D11G12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymidine phosphorylase
Isotype: IgG1	GenBank: BC018160	Gene ID: 1890
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28631	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51270-1 targets PD-ECGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PD-ECGF Monoclonal antibody, PBS Only (Capture) 67104-2-PBS (1B3E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PD-ECGF Monoclonal antibody, PBS Only (Detector) 67104-3-PBS (1D11G12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51270-1, PD-ECGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67104-2-PBS. Detection antibody: 67104-3-PBS. Standard:Ag28631. Range: 1.563-100 ng/mL.