

PD-ECGF Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51270-3**

Capture Antibody Information

Catalog Number: 67104-4-PBS	Clone ID: 1F6D1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymidine phosphorylase
Isotype: IgG1	GenBank: BC018160	Gene ID: 1890
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28631	

Detection Antibody Information

Catalog Number: 67104-6-PBS	Clone ID: 1E6F2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymidine phosphorylase
Isotype: IgG1	GenBank: BC018160	Gene ID: 1890
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28631	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	---------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP51270-3 targets PD-ECGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PD-ECGF Monoclonal antibody, PBS Only (Capture) 67104-4-PBS (1F6D1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PD-ECGF Monoclonal antibody, PBS Only (Detector) 67104-6-PBS (1E6F2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51270-3, PD-ECGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67104-4-PBS. Detection antibody: 67104-6-PBS. Standard:Ag28631. Range: 1.563-100 ng/mL.