

IL22RA1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51277-3**

Capture Antibody Information

Catalog Number: 60866-2-PBS	Clone ID: 1F8H3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 22 receptor, alpha 1
Isotype: IgG1	GenBank: BC029273	Gene ID: 58985
Purification Method: Protein G purification	Immunogen Catalog Number: Ag29983	

Detection Antibody Information

Catalog Number: 60866-4-PBS	Clone ID: 2E12B4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 22 receptor, alpha 1
Isotype: IgG2b	GenBank: BC029273	Gene ID: 58985
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag29983	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51277-3 targets IL22RA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL22RA1 Monoclonal antibody, PBS Only (Capture/Detector) 60866-2-PBS (1F8H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL22RA1 Monoclonal antibody, PBS Only (Detector) 60866-4-PBS (2E12B4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51277-3, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-2-PBS. Detection antibody: 60866-4-PBS. Standard:Ag29983. Range: 3.125-200 ng/mL.