

Moesin Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51280-4**

Capture Antibody Information

Catalog Number:
60868-3-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
1A7F8
Reactivity:
human
GenBank:
BC017293
Immunogen Catalog Number:
Ag23543

Conjugate:
Unconjugated
Full name:
moesin
Gene ID:
4478

Detection Antibody Information

Catalog Number:
60868-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B5E8
Reactivity:
human
GenBank:
BC017293
Immunogen Catalog Number:
Ag23543

Conjugate:
Unconjugated
Full name:
moesin
Gene ID:
4478

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51280-4 targets Moesin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Moesin Monoclonal antibody, PBS Only (Capture) 60868-3-PBS (1A7F8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Moesin Monoclonal antibody, PBS Only (Detector) 60868-5-PBS (1B5E8). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51280-4, Moesin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60868-3-PBS. Detection antibody: 60868-5-PBS. Standard:Ag23543. Range: 1.563-100 ng/mL.