

DLX1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51286-3**

Capture Antibody Information

Catalog Number: 60872-4-PBS	Clone ID: 1H3B4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: distal-less homeobox 1
Isotype: IgG1	GenBank: BC036189	Gene ID: 1745
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28661	

Detection Antibody Information

Catalog Number: 60872-3-PBS	Clone ID: 2F1E11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: distal-less homeobox 1
Isotype: IgG1	GenBank: BC036189	Gene ID: 1745
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28661	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51286-3 targets DLX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DLX1 Monoclonal antibody, PBS Only (Capture) 60872-4-PBS (1H3B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: DLX1 Monoclonal antibody, PBS Only (Detector) 60872-3-PBS (2F1E11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51286-3, DLX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60872-4-PBS. Detection antibody: 60872-3-PBS. Standard:Ag28661. Range: 0.195-100 ng/mL.