

SLC10A1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51290-2**

Capture Antibody Information

Catalog Number:
60876-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1A8C10
Reactivity:
human
GenBank:
BC074724
Immunogen Catalog Number:
Ag25680

Conjugate:
Unconjugated
Full name:
solute carrier family 10 (sodium/bile acid cotransporter family), member 1
Gene ID:
6554

Detection Antibody Information

Catalog Number:
60876-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1A7B7
Reactivity:
human
GenBank:
BC074724
Immunogen Catalog Number:
Ag25680

Conjugate:
Unconjugated
Full name:
solute carrier family 10 (sodium/bile acid cotransporter family), member 1
Gene ID:
6554

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51290-2 targets SLC10A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC10A1 Monoclonal antibody, PBS Only (Capture) 60876-1-PBS (1A8C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC10A1 Monoclonal antibody, PBS Only (Detector) 60876-3-PBS (1A7B7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51290-2, SLC10A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60876-1-PBS. Detection antibody: 60876-3-PBS. Standard:Ag25680. Range: 0.195-100 ng/mL.