

LILRB5 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51294-1**

Capture Antibody Information

Catalog Number: 60879-1-PBS	Clone ID: 1G7D9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5
Isotype: IgG1	GenBank: BC025704	Gene ID: 10990
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28624	

Detection Antibody Information

Catalog Number: 60879-2-PBS	Clone ID: 1F10G10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5
Isotype: IgG1	GenBank: BC025704	Gene ID: 10990
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28624	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51294-1 targets LILRB5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LILRB5 Monoclonal antibody, PBS Only (Capture) 60879-1-PBS (1G7D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: LILRB5 Monoclonal antibody, PBS Only (Detector) 60879-2-PBS (1F10G10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51294-1, LILRB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60879-1-PBS. Detection antibody: 60879-2-PBS. Standard:Ag28624. Range: 0.391-100 ng/mL.