

LILRB5 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51294-2**

Capture Antibody Information

Catalog Number:
60879-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G8F12
Reactivity:
human
GenBank:
BC025704
Immunogen Catalog Number:
Ag28624

Conjugate:
Unconjugated
Full name:
leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5
Gene ID:
10990

Detection Antibody Information

Catalog Number:
60879-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F10G10
Reactivity:
human
GenBank:
BC025704
Immunogen Catalog Number:
Ag28624

Conjugate:
Unconjugated
Full name:
leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5
Gene ID:
10990

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51294-2 targets LILRB5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LILRB5 Monoclonal antibody, PBS Only (Capture) 60879-3-PBS (1G8F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: LILRB5 Monoclonal antibody, PBS Only (Detector) 60879-2-PBS (1F10G10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51294-2, LILRB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60879-3-PBS. Detection antibody: 60879-2-PBS. Standard:Ag28624. Range: 0.391-100 ng/mL.