

PHF17 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51295-2**

Capture Antibody Information

Catalog Number: 60880-3-PBS	Clone ID: 1G7D7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: PHD finger protein 17
Isotype: IgG1	GenBank: BC032376	Gene ID: 79960
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29699	

Detection Antibody Information

Catalog Number: 60880-4-PBS	Clone ID: 1A1A6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: PHD finger protein 17
Isotype: IgG1	GenBank: BC032376	Gene ID: 79960
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29699	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51295-2 targets PHF17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PHF17 Monoclonal antibody, PBS Only (Capture) 60880-3-PBS (1G7D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PHF17 Monoclonal antibody, PBS Only (Detector) 60880-4-PBS (1A1A6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51295-2, PHF17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60880-3-PBS. Detection antibody: 60880-4-PBS. Standard:Ag29699. Range: 0.781-100 ng/mL.