

SPATS2L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51296-1**

Capture Antibody Information

Catalog Number: 60881-1-PBS	Clone ID: 1C5F4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: viral DNA polymerase-transactivated protein 6
Isotype: IgG1	GenBank: BC018736	Gene ID: 26010
Purification Method: Protein G purification	Immunogen Catalog Number: Ag9647	

Detection Antibody Information

Catalog Number: 60881-2-PBS	Clone ID: 1E7A8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: viral DNA polymerase-transactivated protein 6
Isotype: IgG1	GenBank: BC018736	Gene ID: 26010
Purification Method: Protein G purification	Immunogen Catalog Number: Ag9647	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51296-1 targets SPATS2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPATS2L Monoclonal antibody, PBS Only (Capture) 60881-1-PBS (1C5F4). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPATS2L Monoclonal antibody, PBS Only (Detector) 60881-2-PBS (1E7A8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51296-1, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-2-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL.