

GREM1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51302-1**

Capture Antibody Information

Catalog Number: 60888-1-PBS	Clone ID: 1B2C7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: gremlin 1, cysteine knot superfamily, homolog (Xenopus laevis)
Isotype: IgG1	GenBank: BC093778	Gene ID: 26585
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag12734	

Detection Antibody Information

Catalog Number: 60888-2-PBS	Clone ID: 1B9B6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: gremlin 1, cysteine knot superfamily, homolog (Xenopus laevis)
Isotype: IgG1	GenBank: BC093778	Gene ID: 26585
Purification Method: Protein G purification	Immunogen Catalog Number: Ag12734	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51302-1 targets GREM1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GREM1 Monoclonal antibody, PBS Only (Capture) 60888-1-PBS (1B2C7). 100 µg. Concentration 1 mg/mL.

Detection antibody: GREM1 Monoclonal antibody, PBS Only (Detector) 60888-2-PBS (1B9B6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51302-1, GREM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60888-1-PBS. Detection antibody: 60888-2-PBS. Standard:Ag12734. Range: 0.098-100 ng/mL.