

GOLM1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51312-1**

Capture Antibody Information

Catalog Number: 66331-2-PBS	Clone ID: 1G8E7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: golgi membrane protein 1
Isotype: IgG1	GenBank: BC001740	Gene ID: 51280
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7207	

Detection Antibody Information

Catalog Number: 66331-3-PBS	Clone ID: 2G6H3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: golgi membrane protein 1
Isotype: IgG1	GenBank: BC001740	Gene ID: 51280
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7207	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51312-1 targets GOLM1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GOLM1 Monoclonal antibody, PBS Only (Capture) 66331-2-PBS (1G8E7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: GOLM1 Monoclonal antibody, PBS Only (Capture/Detector) 66331-3-PBS (2G6H3). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51312-1, GOLF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66331-2-PBS. Detection antibody: 66331-3-PBS. Standard:Ag7207. Range: 0.098-100 ng/mL.