

ITLN1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51319-1**

Capture Antibody Information

Catalog Number: 60906-1-PBS	Clone ID: 1E9H1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: intelectin 1 (galactofuranose binding)
Isotype: IgG3	GenBank: BC020664	Gene ID: 55600
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag28617	

Detection Antibody Information

Catalog Number: 60906-2-PBS	Clone ID: 3A3D7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: intelectin 1 (galactofuranose binding)
Isotype: IgG1	GenBank: BC020664	Gene ID: 55600
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28617	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51319-1 targets ITLN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ITLN1 Monoclonal antibody, PBS Only (Capture) 60906-1-PBS (1E9H1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ITLN1 Monoclonal antibody, PBS Only (Detector) 60906-2-PBS (3A3D7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS.



Cytometric bead array standard curve of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS. Standard:Ag28617. Range: 0.098-100 ng/mL