

HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51320-1**

Capture Antibody Information

Catalog Number: 67445-2-PBS	Clone ID: 3D2D10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: heterogeneous nuclear ribonucleoprotein A2/B1
Isotype: IgG1	GenBank: BC000506	Gene ID: 3181
Purification Method: Protein G purification	Immunogen Catalog Number: Ag6699	

Detection Antibody Information

Catalog Number: 67445-3-PBS	Clone ID: 1C2B7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: heterogeneous nuclear ribonucleoprotein A2/B1
Isotype: IgG1	GenBank: BC000506	Gene ID: 3181
Purification Method: Protein G purification	Immunogen Catalog Number: Ag6699	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51320-1 targets HNRNPA2B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HNRNPA2B1 Monoclonal antibody, PBS Only (Capture) 67445-2-PBS (3D2D10). 100 µg. Concentration 1 mg/mL.

Detection antibody: HNRNPA2B1 Monoclonal antibody, PBS Only (Detector) 67445-3-PBS (1C2B7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51320-1, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-3-PBS.



Cytometric bead array standard curve of MP51320-1, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-3-PBS. Standard:Ag6699. Range: 0.098-100 ng/mL.