

PRSS8 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51333-1**

Capture Antibody Information

Catalog Number: 60912-1-PBS	Clone ID: 2E8G7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: protease, serine, 8
Isotype: IgG1	GenBank: BC001462	Gene ID: 5652
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7234	

Detection Antibody Information

Catalog Number: 60912-2-PBS	Clone ID: 1H2A7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: protease, serine, 8
Isotype: IgG1	GenBank: BC001462	Gene ID: 5652
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7234	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51333-1 targets PRSS8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRSS8 Monoclonal antibody, PBS Only (Capture) 60912-1-PBS (2E8G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PRSS8 Monoclonal antibody, PBS Only (Detector) 60912-2-PBS (1H2A7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51333-1, PRSS8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60912-1-PBS. Detection antibody: 60912-2-PBS.



Cytometric bead array standard curve of MP51333-1, PRSS8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60912-1-PBS. Detection antibody: 60912-2-PBS. Standard:Ag7234. Range: 3.125-200 ng/mL