

IGSF21 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51337-2**

Capture Antibody Information

Catalog Number: 60915-1-PBS	Clone ID: 2A8C5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: immunoglobulin superfamily, member 21
Isotype: IgG2a	GenBank: BC007618	Gene ID: 84966
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag15062	

Detection Antibody Information

Catalog Number: 60915-3-PBS	Clone ID: 1F4E6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: immunoglobulin superfamily, member 21
Isotype: IgG1	GenBank: BC007618	Gene ID: 84966
Purification Method: Protein G purification	Immunogen Catalog Number: Ag15062	

Applications

Tested Applications: Cytometric bead array	Range: 1.953-500 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51337-2 targets IGSF21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IGSF21 Monoclonal antibody, PBS Only (Capture) 60915-1-PBS (2A8C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: IGSF21 Monoclonal antibody, PBS Only (Detector) 60915-3-PBS (1F4E6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51337-2, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-3-PBS.



Cytometric bead array standard curve of MP51337-2, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-3-PBS. Standard: Ag15062. Range: 1.953-500 ng/mL.