

SFRS2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51340-1**

Capture Antibody Information

Catalog Number:
60922-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2A8G8
Reactivity:
human
GenBank:
BC070086
Immunogen Catalog Number:
Ag17177

Conjugate:
Unconjugated
Full name:
splicing factor, arginine/serine-rich 2
Gene ID:
6427

Detection Antibody Information

Catalog Number:
60922-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C5A9
Reactivity:
human
GenBank:
BC070086
Immunogen Catalog Number:
Ag17177

Conjugate:
Unconjugated
Full name:
splicing factor, arginine/serine-rich 2
Gene ID:
6427

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51340-1 targets SFRS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SFRS2 Monoclonal antibody, PBS Only (Capture) 60922-1-PBS (2A8G8). 100 µg. Concentration 1 mg/mL.

Detection antibody: SFRS2 Monoclonal antibody, PBS Only (Detector) 60922-2-PBS (1C5A9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51340-1, SFRS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60922-1-PBS. Detection antibody: 60922-2-PBS.



Cytometric bead array standard curve of MP51340-1, SFRS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60922-1-PBS. Detection antibody: 60922-2-PBS. Standard: Ag17177. Range: 3.125-400 ng/mL.