

# FOXO3A Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51341-1**

## Capture Antibody Information

**Catalog Number:**  
66428-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1F12D11  
**Reactivity:**  
human, mouse, rat  
**GenBank:**  
BC021224  
**Immunogen Catalog Number:**  
Ag1289

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box O3  
**Gene ID:**  
2309

## Detection Antibody Information

**Catalog Number:**  
66428-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
1D11G4  
**Reactivity:**  
human  
**GenBank:**  
BC021224  
**Immunogen Catalog Number:**  
Ag1289

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box O3  
**Gene ID:**  
2309

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.391-400 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51341-1 targets FOXO3A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXO3A Monoclonal antibody, PBS Only (Capture) 66428-1-PBS (1F12D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: FOXO3A Monoclonal antibody, PBS Only (Detector) 66428-2-PBS (1D11G4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

Selected Validation Data



Sample test of MP51341-1, FOXO3A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66428-1-PBS. Detection antibody: 66428-2-PBS.



Cytometric bead array standard curve of MP51341-1, FOXO3A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66428-1-PBS. Detection antibody: 66428-2-PBS. Standard:Ag1289. Range: 0.391-400 ng/mL