

KIAA0494 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51343-1**

Capture Antibody Information

Catalog Number: 60924-1-PBS	Clone ID: 1F5E5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: KIAA0494
Isotype: IgG2a	GenBank: BC002525	Gene ID: 9813
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag33197	

Detection Antibody Information

Catalog Number: 60924-2-PBS	Clone ID: 1H1B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: KIAA0494
Isotype: IgG1	GenBank: BC002525	Gene ID: 9813
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33197	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51343-1 targets KIAA0494 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA0494 Monoclonal antibody, PBS Only (Capture) 60924-1-PBS (1F5E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: KIAA0494 Monoclonal antibody, PBS Only (Detector) 60924-2-PBS (1H1B5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51343-1, KIAA0494 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60924-1-PBS. Detection antibody: 60924-2-PBS.



Cytometric bead array standard curve of MP51343-1, KIAA0494 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60924-1-PBS. Detection antibody: 60924-2-PBS. Standard:Ag33197. Range: 0.391-400 ng/mL.