

RBM3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51344-2**

Capture Antibody Information

Catalog Number:
60925-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D2G1

Reactivity:
human

GenBank:
BC006825

Immunogen Catalog Number:
Ag28329

Conjugate:
Unconjugated

Full name:
RNA binding motif (RNP1, RRM)
protein 3

Gene ID:
5935

Detection Antibody Information

Catalog Number:
60925-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E8E1

Reactivity:
human

GenBank:
BC006825

Immunogen Catalog Number:
Ag28329

Conjugate:
Unconjugated

Full name:
RNA binding motif (RNP1, RRM)
protein 3

Gene ID:
5935

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51344-2 targets RBM3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBM3 Monoclonal antibody, PBS Only (Capture) 60925-1-PBS (2D2G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBM3 Monoclonal antibody, PBS Only (Detector) 60925-3-PBS (1E8E1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51344-2, RBM3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60925-1-PBS. Detection antibody: 60925-3-PBS.



Cytometric bead array standard curve of MP51344-2, RBM3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60925-1-PBS. Detection antibody: 60925-3-PBS. Standard:Ag28329. Range: 0.391-50 ng/mL.