

Aiolos/IKZF3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51356-1**

Capture Antibody Information

Catalog Number:
60938-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3D7E2
Reactivity:
human
GenBank:
BC032707
Immunogen Catalog Number:
Ag5537

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 3 (Aiolos)
Gene ID:
22806

Detection Antibody Information

Catalog Number:
60938-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
1G8A4
Reactivity:
human
GenBank:
BC032707
Immunogen Catalog Number:
Ag5537

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 3 (Aiolos)
Gene ID:
22806

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51356-1 targets Aiolos/IKZF3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Aiolos/IKZF3 Monoclonal antibody, PBS Only (Capture) 60938-1-PBS (3D7E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Aiolos/IKZF3 Monoclonal antibody, PBS Only (Detector) 60938-2-PBS (1G8A4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51356-1, Aiolos/IKZF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60938-1-PBS. Detection antibody: 60938-2-PBS.



Cytometric bead array standard curve of MP51356-1, Aiolos/IKZF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60938-1-PBS. Detection antibody: 60938-2-PBS. Standard: Ag5537. Range: 0.391-400 ng/mL.